

## Acetone precipitation of proteins

---

### Reagents & Solutions

- Ice-cold Acetone, HPLC grade (-20°C) \*

### Equipment

- Benchtop Eppendorf centrifuge \*

\* available in the core facility

### Protocol

1. Add 4x volume of ice-cold acetone to the protein sample<sup>1</sup>
2. Incubate for 15 min at -80°C, followed by 90-120 min incubation at -20°C<sup>2</sup>
3. Centrifuge the tubes for 15 min at 16.000 × g
4. Carefully remove the supernatant<sup>3</sup>
5. Add 500 µL ice-cold acetone to the protein pellet
6. Centrifuge for 5 min at 16.000 × g, discard supernatant
7. Repeat this step one more time
8. Air-dry pellet and resuspend in buffer of choice, depending on the digestion protocol

## Notes

<sup>1</sup> Use 2 mL tubes for sample volumes larger than 250  $\mu$ L. Use 5 mL tubes for sample volumes larger than 350  $\mu$ L. Use 10 mL falcon tubes for samples larger than 900  $\mu$ L.

<sup>2</sup> Precipitation overnight is also possible, but not necessary.

<sup>3</sup> In case of low amounts of proteins it can happen that the pellet is not clearly visible. Add GlycoBlue co-precipitant to increase visibility of the proteins.